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
A handwritten signature in blue ink, consisting of a series of loops and curves, positioned above the name Dr. M.J. van de Velde.

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Date: September 1, 2020
Subject: Registration In-vitro diagnostics

Dear Mr. Wei,,

On August 28, 2020, I received your notification pursuant to Article 4, paragraph 1 of the Dutch Decree in vitro diagnostics (BIVD) under the company name Anbio (Xiamen) Biotechnology Co., Ltd. with European authorized Lotus NL B.V. market the product below as an in-vitro diagnostic agent product on the European market.

The product is registered as an in vitro diagnostic under number:

COVID-19 Test Kit(Real-time PCR),
Rapid COVID-19 Antigen Test (Colloidal Gold)
(no brand name) (NL-CA002-2020-53206)

With this you have fulfilled your obligation under Article 4 BIVD.

In all further correspondence regarding the above-mentioned product, I request that you mention this number. No further rights can be derived from this number, it only serves to facilitate the administration of the notification.

The registration of in vitro diagnostics as a medical device on the basis of the Classification criteria (Annex II) to Directive 98/79/EC on medical devices for in vitro diagnostics is subject to possible revisions of European regulations on the classification of medical devices and progressive scientific understanding (see Article 10, paragraph 1 of Directive 98/79 / EC).

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Attachments

Your request
28 augustus 2020

Correspondence from close and to the return address, stating the date and mark of this letter.


Notification of in vitro diagnostic medical devices implies that the manufacturer, Anbio (Xiamen) Biotechnology Co., Ltd. has affixed the CE conformity marking to the product concerned before placing it on the market in an EU Member State. In this way, Lotus NL B.V. guarantees that the in vitro diagnostic icum complies with the essential requirements as included in Annex I to Directive 98/79 / EC (and in the corresponding part 1 of the Decree)

For the sake of completeness, we would like to point out that an in-vitro diagnostic must meet the requirements of the BIVD. The BIVD is based on the guideline for in vitro diagnostics, 98/79/EC. In particular, we point out the Dutch language requirement as it applies in the Netherlands, the requirements for keeping the technical documentation available and the obligation to have a Post Marketing Surveillance and vigilanciesysteem.

Finally, I note that your notification - the administrative notification as a manufacturer - and this letter does not constitute a judgment on the status or qualification of your product: notification does not mean that there is actually an in-vitro diagnostic in the sense of the present legislation and regulations. Where appropriate, the Healthcare and Youth Inspectorate (IGJ), charged with supervising compliance with the provisions laid down by or pursuant to the law, may take a position on the status of a product, in accordance with established case law. It is up to the national court to determine whether a product falls within the definition of in vitro diagnostics.

The Minister for Medical Care and Sport,
on their behalf,

Department head
Farmatec

A handwritten signature in blue ink, appearing to be 'M.J. van de Velde', written over a horizontal line.

Dhr. M.J. van de Velde

DECLARATION OF CONFORMITY

According Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex III.

Manufacturer: Anbio (Xiamen) Biotechnology Co.,Ltd.

Address: No.2016, Wengjiao West Road, Xinyang Street, Haicang District,361026 Xiamen, Fujian, China.

European Representative: Lotus NL B.V.

Contact person: Peter **E-mail:** peter@lotusnl.com

Address: Koningin Julianaplein 10,1e Verd, 2595AA,The Hague, Netherlands.

In Vitro Diagnostic Directive:

- **Rapid COVID-19 Antigen Test (Colloidal Gold)**

Category: Others.

Conformity assessment route: Declaration of Conformity IVDD Annex III

Applicable Standards:

ISO 13485:2016

ISO 14971:2019

EN ISO 18113-1:2011

EN ISO 18113-2:2011

EN ISO 18113-3:2011

EN 13641:2002

ISO 15223-1:2016

EN 13612:2002

ISO 23640:2015

EN 62366-1:2015

We, the manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In Vitro Diagnostic Medical Devices.

We agree to develop, implement and maintain a documented post-production monitoring process.

Signed on:04/09/2020

Place:Xiamen,China

Name of authorized signatory: *Daming Wang*

Position held in the company: General Manager

Seal/Stamp:

Anbio (Xiamen) Biotechnology Co.,Ltd.

Anbio Biotech

Rapid COVID-19 Antigen test
(Colloidal Gold)

Anbio (Xiamen) Biotechnology Co.,Ltd.

Product information

- ❖ Registration certificate: CE: NL-CA002-2020-53206
- ❖ Specification: 20tests/kit;
- ❖ Specimen type: Human nasal swabs , throat swabs or saliva
- ❖ Storage and shelf life: 2-30°C ; 24 months
- ❖ Reagent Kit including:


Contents	20 tests/kit
COVID-19 Antigen Test Card	20
Extraction Solution	20
Swab	20
Instruction for Use	1

Detection Principle

This reagent is based on colloidal gold immunochromatography assay.

During the test, specimen extracts are applied to the Test Cards. If there were SARS-CoV-2 antigen in the extract, the antigen will bind to the SARS-CoV-2 monoclonal antibody. During lateral flow, the complex will move along the nitrocellulose membrane toward the end of the absorbent paper. When passing the test line (line T, coated with another SARS-CoV-2 monoclonal antibody) the complex is captured by SARS-CoV-2 antibody on test line resulting in coloring on line T; when passing the line C, colloidal gold-labeled goat anti-rabbit IgG is captured by control line (line C, coated with rabbit IgG) resulting in coloring on line C.

Performance Evaluation——Sensitive result: 98.48 %; Specificity result: 100%

Clinical performance of Rapid COVID-19 Antigen Test(Colloidal Gold) was determined by testing 591 positive and 547 negative specimens

		PCR result		
		Positive	Negative	Total
Rapid COVID-19	Positive	582(a)	0(b)	582(a+b)
Antigen Test(Colloidal	Negative	9(c)	547(d)	556(c+d)
Gold) result	Total	591(a+c)	547(b+d)	1138(a+b+c+d)

Coincidence rate and 95% confidence interval		
	Coincidence rate	95% confidence interval
Clinical sensitivity	98.48%	97.13%~99.30%
Clinical specificity	100%	99.33%~100%
Total coincidence rate	99.21%	98.50%~99.64%

Performance Evaluation——Ability of anti-jamming: 100%

Antigen positive samples for other respiratory virus*13,Extraction solution*10ml(50 test cards)

No	Sample classification	Test result
N1	2019-nCoV virus nucleic acid negative FluA positive	Negative
N2	2019-nCoV virus nucleic acid negative Metapneumo virus positive	Negative
N3	2019-nCoV virus nucleic acid negative FluB positive	Negative
N4	2019-nCoV virus nucleic acid negative parainfluenza virus positive	Negative
N5	2019-nCoV virus nucleic acid negative respiratory syncytial virus positive	Negative
N6	2019-nCoV virus nucleic acid negative rhino virus positive	Negative
N7	2019-nCoV virus nucleic acid negative adeno virus positive	Negative
N8	2019-nCoV virus nucleic acid negative mycoplasma pneumoniae positive	Negative
N9	2019-nCoV virus nucleic acid negative chlamydia pneumoniae positive	Negative
N10	2019-nCoV virus nucleic acid negative Human coronavirus 229E positive Human	Negative
N11	2019-nCoV virus nucleic acid negative Human corona virus OC43 positive	Negative
N12	2019-nCoV virus nucleic acid negative corona virus NL63 positive	Negative
N13	2019-nCoV virus nucleic acid negative Human corona virus HKU1 positive	Negative
Diluent (N14-N64)	/	Negative
Coincidence Rate	100%	

Features and Benefits


***It can be tested without professional lab**

——**Can be used anywhere**


——**No need the professional person**


Three specimen types:

Throat swabs

Nasal swabs

Saliva

Features and Benefits

***No equipment required**

- The result can be read by eyes
- Easy to use

Features and Benefits

*Fast test result:

——15mins

——No need waiting for a long time

——Screening a large number of persons in a short time


Features and Benefits

*Stable storage: At 2-30°C for 24 months

—— Easy to shipping and storage


—— No restrictions on transport, Can be shipping by Plane/ Van/Train/Ship, etc

—— No need to refrigerate

Specimen Collection and Preparation


Throat Swab Specimen Collection:


Let the patient's head tilt slightly, mouth open, and make "ah" sounds, exposing the pharyngeal tonsils on both sides. Hold the swab and wipe the pharyngeal tonsils on both sides of the patient with moderate force back and forth for at least 3 times.

Nasal Swab Specimen Collection:


1. Insert the swab into one nostril of the patient. The swab tip should be inserted up to 2.5 cm (1 inch) from the edge of the nostril.
2. Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected
3. Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities. Withdraw the swab from the nasal cavity.

Cough deeply, make the noise of "kruuuu"


Place the swab on your tongue

Hold the swab in your mouth for 10-20s, moisten the swab with saliva


Test Method

The test should be operated at room temperature (15-30°C) .

1. Place the extraction tube with opening facing up.
2. Put the swab with collected specimen into the extraction tube, hold and press the swab head against the wall of tube with force while rotating the swab for about 10 seconds (10-15 times) to release the antigen into the extraction solution from the swab head.
3. Removing swab: Squeeze the swab head while removing the swab in order to remove as much liquid as possible from the swab. Dispose of swabs according to biohazard waste disposal regulations.
4. Screw the nozzle cap onto the extraction tube.
5. Loading: drip 2 drops of swab specimen or 3 drops of saliva specimen into the sample well of the Test Card, and start the timer.
6. Read the results at 15minutes.

Interpretation of test results

【Interpretation of test results】

Negative (-):

There is coloration on line C only showing as following picture, suggesting that there is no SARS-CoV-2 antigen in the specimen.

Positive (+):

There are coloration on both line C and line T showing as follow pictures, suggesting that there is SARS-CoV-2 antigen in the specimen.

Invalid:

There is no coloration on line C, as shown in the following pictures. The test is invalid or an error in operation occurred. Repeat the assay with a new cartridge

About ANBIO

*Anbio Biotechnology Co.,Ltd. founded in 2015 , is located in Xiamen, China.

*Specialized in rapid test and production of POCT in vitro diagnostic devices and reagents

*Anbio takes POCT as its main development direction and focuses on the accuracy, rapidity and intelligence of clinical diagnosis.

*Products cover immune, biochemical, blood grouping, molecular diagnosis and other technical platforms, and through product integration, the company is fully More than 70% of medical laboratory needs.

About ANBIO

Anbio has a senior R & D, production, management and operation team; a full industrial chain development layout from biological raw materials, medical laboratory products to professional services, Project cooperation with Suzhou Institute of Medical Engineering, Chinese Academy of Sciences, School of Environment and Public Health, Huasha College and other institutions

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
Dr. M.J. van de Velde

Pagina 2 van 2

ISO 13485 and Medical Equipment Manufacturers Licence

BUREAU VERITAS
Certification

Anbio (Xiamen) Biotechnology Co., Ltd
No.2016, Wengjiao West Road, Xinyang Street, Haicang District, Xiamen City,
P.R. China

Certified site:
NO.2016, WENGJIAO WEST ROAD, XINYANG STREET, HAICANG DISTRICT, XIAMEN CITY,
P.R. CHINA

*Bureau Veritas Italia S.p.A. certifies that the Management System of the
above organisation has been audited and found to be in accordance
with the requirements of the management system standards detailed below*

EN ISO 13485:2016
Scope of certification

Design and manufacture of fluorescence-based immunoassay reagent
and devices, as an aid in clinical assessment of cardiovascular,
gastric, inflammation, diabetic and infectious diseases detection, as
well as hormone, vitamin testing. Design and manufacture of IVD
reagents as an aid in clinical assessment of blood type testing.

Certificate awarded in conformity with the requirements of ACCREDIA DT 02-DC Rev.00

Original cycle start date:	22/06/2020
Expiry date of previous cycle:	n.a.
Certification / Recertification Audit date:	17/05/2020
Certification / Recertification cycle start date:	22/06/2020

Subject to the continued satisfactory operation of the organization's
Management System, this certificate expires on: **21/06/2023**

Certificate No. - Version: **IT298645-1** Revision date: **22/06/2020**

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Further clarifications regarding the scope of this certificate and the applicability of the
management system requirements may be obtained by consulting the organisation.
To check this certificate validity please refer to the website www.bureauveritas.it

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Packing information of Rapid COVID-19 Antigen test

1 Package for Test Kit

Item	Package for Test Kit
Image	
Size (L*W*H,CM)	18*14*7.5 cm
Weight	310 g
Specifications	20 Tests/Kit

2 Package for Carton

Item	Package for Test Kit
Image	
Size (L*W*H,CM)	72*38*40cm
Weight	18Kg
Volumetric weight	≈20kg
Specifications	1000 Tests (50 Kits) /Carton

3 Main components

No.	Components	Image
1	Test Card*20	
2	Extraction solution*20	
3	Swab*20	
4	Instruction for Use*1	

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Bureau Veritas Italia S.p.A., Viale Monza, 347 - 20126 Milano, Italia